

Technical data sheet

Attachment no. 6 pt. 1 subpt. 1.1

Edition II 2019.02.27

Page 1/2

TD-41 dicoNEX Single use syringe, 3-part

| MEDICAL DEVICE DATA: | | | |
|---|---|--|-------|
| CATALOGUE NUMBER | 010ML-3CZ-BL 010ML-3CZ-LL-BL | | |
| MANUFACTURER | ZARYS International Group sp. z o.o. sp.k., Poland | | |
| CLASSIFICATION | IIa | RULE | 2 |
| METHOD OF STERILIZATION | Ethylene oxide | CODE OF GENERIC DEVICE GROUP ACCORDING TO GMDN | 47017 |
| PERFORMANCE CHARACTERISTICS | 3-part syringe for intramuscular, intravenous and subcutaneous injection Intended for administering drugs / liquids / food to the patient. It can also be used to take body fluid samples from the patient | | |
| CURRENT TECHNICAL STANDARDS ARE APPLICABLE TO THE DEVICE | ISO 7886-1 Sterile hypodermic syringes for single use | | |

| TECHNICAL DATA: | | | |
|-----------------|--|--|--|
| TECHNICAL DATA. | | | |
| RAW MATERIALS | barrel – polypropylene plunger/piston – polypropylene rubber sealing – polyisoprene rubber (synthetic) silicone oil (polysiloxane) | | |
| COLOUR | barrel - transparent plunger/piston - green CAL 19.01.2021 rubber sealing - black silicone oil - transparent | | |
| FEATURES | 4. silicone oil – transparent three-piece syringe transparent barrel allows perfect visualization of the content green contrasting plunger luer or luer-lock nozzle with dead space 0,035 ml sealing in the form of a double ring on a plunger stopper situated on the top of plunger smooth plunger movement thanks to rubber sealing and silicone oil on the inner surface of the barrel limitation flange securing against accidental removal of the plunger name of the syringe and manufacturer's logo placed on the barrel clear, legible and permanent black graduation facilitating dosing ever 0,2ml single-use, non-pyrogenic, non-toxic latex-free phthalate-free EO sterilized shelf life: 5 years | | |
| VARIABLES | 010ML-3CZ-BL with luer nozzle 010ML-3CZ-LL-BL with luer-lock nozzle | | |
| SIZE/CAPACITY | 10 ml | | |



Technical data sheet

Attachment no. 6 pt. 1 subpt. 1.1

Edition II 2019.02.27

Page 2/2

TD-41 dicoNEX Single use syringe, 3-part

