

## Fiche technique : Seringues 2 pièces DuoNEX

### duoNEX

single-use syringe, 2-part, Luer

#### features

- white line
- transparent barrel allows full visual control of content
- white contrasting plunger
- Luer nozzle centrally located in 2 ml syringe, eccentric in 5, 10, 20 ml syringes
- smooth plunger sliding
- limitation flange securing against accidental removal of plunger
- clear, legible and permanent black graduation facilitating dosage, adequately:
  - every 0,1 ml for capacity of 2 ml
  - every 0,2 ml for capacity of 5 ml
  - every 0,5 ml for capacity of 10 ml
  - every 1 ml for capacity of 20 ml
- extended graduation that includes a 10% extension of the nominal capacity
- name of syringe and manufacturer's logo placed on barrel
- non-pyrogenic, non-toxic
- packaging: 1 pc./paper-foil

#### indication

- 2-piece syringe for subcutaneous, intramuscular and intravenous infusions intended for drawing up and administering of medicines/liquids to patient and collecting body fluid samples from patient



| REF         | type | capacity | position of nozzle | intermediate packaging | bulk packaging |
|-------------|------|----------|--------------------|------------------------|----------------|
| 002ML-2CZ-B | Luer | 2 ml     | centric            | 100 pcs                | 30 x 100 pcs   |
| 005ML-2CZ-B | Luer | 5 ml     | eccentric          | 100 pcs                | 24 x 100 pcs   |
| 010ML-2CZ-B | Luer | 10 ml    | eccentric          | 100 pcs                | 16 x 100 pcs   |
| 020ML-2CZ-B | Luer | 20 ml    | eccentric          | 50 pcs                 | 16 x 50 pcs    |

